



# Recombinant *Xenopus tropicalis* MAP kinase-interacting serine/threonine-protein kinase 2 (mknk2)

<b>Product Code</b>	CSB-MP724126XBF
<b>Abbreviation</b>	mknk2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q66I46
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVQKKTTEMK GFHRSFKGQN PFDAAYEMES RNMASVFNFD CPSRPDVPSS APIDIPDAKK RTKKKKRCRA TDSFTGRFED MYQLQQEILG EGAYAKVQSC INLITNKEYA VKIIEKRP GH SRSRVFREVE MLYQCQGH SN VLELIEFFEE EDKFYLVFEK MCGGSILNHI HRRRHFNERE ASFVVRDIAE ALNYLHNKGI AHRDLKPENI LCESPHQVSP VKICDFDLGS GIKLNSDCSP ISTPELLTPC GSAEYMAPEV VEA FN EEASI YDKRCDLWSL GVILYIMLSG YPPFVGHC GS DCGWDRGEAC PACQNMLFVS IQEGKYE FPE KDWAHISSGA KDLISKLLLR DAKKRLSAAQ VLQHPWVQGN APDNTLPTPI ILQRNSSAKD LTSFAAEAIA MNRQLMERE EEEGTESSSS CPIVVKATSC SMQLSP PSES KLA KR RHGSK GGISPPSLAP LLIVSD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	mknk2
<b>Protein Names</b>	Recommended name: MAP kinase-interacting serine/threonine-protein kinase 2 EC= 2.7.11.1 Alternative name(s): MAP kinase signal-integrating kinase 2 Short name= MAPK signal-integrating kinase 2 Short name= Mnk2
<b>Expression Region</b>	1-466
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.