



Recombinant *Xenopus laevis* Charged multivesicular body protein 2b-B (chmp2b-b)

Product Code	CSB-MP724131XBE
Abbreviation	chmp2b-b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66IV6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MASLFKKKTVDIIIREQNKE LRGTQRAITR DRAALEKQEK QLEMEIKKMA KTGNKDACRV LAKQLVQLRK QKTRTYAVSS KVTSMSTQTK VMSSQMKMAG AMSTTAKTMQ AVNKKMDPQK TLQTMQNFQK ENMKMEMTDE MINDTLDDIF DASEDEEESQ DIVNQVLDEI GIEISGKMAK APSAAKGLPS TSAAKSKGIS DEEIERQLKA LGVD
Source	Mammalian cell
Target Names	chmp2b-b
Protein Names	Recommended name: Charged multivesicular body protein 2b-B Alternative name(s): Chromatin-modifying protein 2b-B Short name= CHMP2b-B
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.