



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 69

(Os09g0459600, LOC_Os09g28560)

Product Code	CSB-BP724185OFG
Abbreviation	Os09g0459600, LOC_Os09g28560
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67J17
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAMTAAAVTV PLGVLLRREV TSEMERPDV LCGEAARSRK GEDFTLLLAE AGERVAGDPS TSFSVFALFD GHNGSGAAMY AKKNLLNLL RAIPSGLSRD EWLAVLPRAL VAAFVKTDKD FQAVAETSGT TVTFVVIDEW VVTVASVGDS RCILESADGS LYHLSADHRF DSNQDEVQRV TACGSKVGKL NLVGGPEVGP LRCWPGGLCL SRSIGDMDVG ECIIPVPHVK QVKLSNAGGR IIIASDGVWD DLTFEMALEC SRGFPSDIAA NRIVNEAIHP RGLRDDTTTCI VVDILPPEKL APSPPTKRQG KIVFNNMFRR KHTDVSFILD REYAEPDEVE EIFDDGSAML SKRLAAGYAL QSMFEPFSCA VCQVQLKAGQ GISVHSNPLQ HEKLQGWQGP FLCQSCNEKK DAIEGKRPPR DS
Source	Baculovirus
Target Names	Os09g0459600
Protein Names	Recommended name: Probable protein phosphatase 2C 69 Short name= OsPP2C69 EC= 3.1.3.16
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.