



Recombinant Gamma-glutamyl phosphate reductase (proA)

Product Code	CSB-MP724200SRW
Abbreviation	proA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67LC2
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MDVRSMARQA RAAQRQLELT SPEERNAALE AIARALEARA GQIVAANRAD LAAEEAGLP GPMIARLRD EGKLAGVIRG VRAVAELPDP LAVEPQAWIR PNGLRIQVR VPLGVVGIIY ESRPGVTVD AVLCLKAGNA VLLKGGKEAA RSNAALGEAM AAGLQEVGLP VELAQVLPST REAAQELMAA RSLVDVLIPR GGAGLIRATV EQSQVPVIET GTGVCHVYVH RAADLAMARE IVLNAKCSNP AVCNAAETLL VDREVAAPFL AEAGPALLAA GVRLRACPEA LAILQETARP DLAGAGPAAV EPATEEDWAA EYLDLCLAVK VVSGLDEALA HIARYGTGHS EAIVTGDAGA ADRFLRSVDA AAVYHNASTR FTDGGEFGFG AEIGISTQKL HARGPMGLAE LTTYKYVVLG NGQVR
Source	Mammalian cell
Target Names	proA
Protein Names	Recommended name: Gamma-glutamyl phosphate reductase Short name= GPR EC= 1.2.1.41 Alternative name(s): Glutamate-5-semialdehyde dehydrogenase Glutamyl-gamma-semialdehyde dehydrogenase Short name= GSA dehydrogenase
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.