



Recombinant Tyrosine--tRNA ligase (tyrS)

Product Code	CSB-EP724218SRW
Abbreviation	tyrS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67QD6
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MTMAFMSVDE QMKILMRGVV DLVSEELRQ KLERSVKTGR PLRVKLGIDP TGKDLTLGHT VPLRKLRDFV ECGHQGVLLI GDYTAMVGDP TGRNEARPQL THAETTANAQ RYLEQAARVL DVSKLEIRRN SEWLAPMSFS DVIRLAAKST VARMLEREDF KKRYEEGRPI FIHEFFYPLM QGTDSDVAVQA DVELGGTDQK FNLLAGRDQ RDAGQEPQVC LMTPIVEGLD GVQKMSKSLG NYIGLDHTPD EMFGRMTSIP DSLIITYFTY FTDVPQEEIE RIEAAMAAGE NPMTFKKQLG RAITTYGGT EEEARLAEER WVAQFSRGEV PEDIPDVVLP AAELPMQAAR VLFTAGLAPS LSEARRLIEQ GGFTVDGEKV TDPRAELALR PGQVLKAGKR KYGRVVVK
Source	E.coli
Target Names	tyrS
Protein Names	Recommended name: Tyrosine--tRNA ligase EC= 6.1.1.1 Alternative name(s): Tyrosyl-tRNA synthetase Short name= TyrRS
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.