



Recombinant Danio rerio WD repeat-containing protein 18 (wdr18)

Product Code	CSB-EP724253DIL-B
Abbreviation	wdr18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68E10
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSAPVEVVLS ADSAGQLFNC AVYDPHTGSE FLSYRGGNTS SRSLTILNGE YILGAQLGKN YINVWEIQRK DQLQQKIVCP GIVTCLCASP DGLYVLAGIA EAIYLWEVST GNLLAILSRH FQDLSCIKFT DDSSHFVSGG KDNLAFIWNL SSVVQLDSSR TPEPRHILSR HSLPITDIHC GLMGPQARVA TASLDQTVKV WEISSGEMLL SVLFDVGIMS VTFDPCYFL FCGGSDGNIF QVSLCSTSL RDKTFQSDSE GNQVFKGHRN LVTCLSVSMD GTVLLSGSND ETVRMWDVQS KQCIWTINHR GPVTNAAIIP APANMFLSDS HPAVPLPRFS RHLNPSEQGD GTGTGGMSLR LGANTQEPEG TYLEKAEELY SLMCAVTDKS VFGDGENTKV RVSELEEEVR TLKKINKDLY EFSTQLLTKP N
Source	E.coli
Target Names	wdr18
Protein Names	Recommended name: WD repeat-containing protein 18
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.