



Recombinant Human herpesvirus 6A Protein RJ1 (RJ1)

Product Code	CSB-BP724341HJZ
Abbreviation	RJ1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69577
Product Type	Recombinant Protein
Immunogen Species	Human herpesvirus 6A (strain Uganda-1102) (HHV-6 variant A) (Human B lymphotropic virus)
Purity	≥85% (SDS-PAGE)
Sequence	RFSV RDSFMVSVIV GVRVRVDVLV RVNVLVRARP RVRVRARPRV RVRVRARPRV RVRVRASPRV RVRVITPPFF SPPGGFDVIF KTRGIPSPGE TVLDKKVSRP WRGLTSRRGA TMI
Source	Baculovirus
Target Names	RJ1;
Protein Names	Recommended name: Protein RJ1 Alternative name(s): Protein LT1
Expression Region	27-143
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.