



Recombinant *Oryza sativa* subsp. *japonica* Ribonuclease 3-like protein 3 (Os06g0358800, LOC_Os06g25250)

Product Code	CSB-BP724351OFG
Abbreviation	Os06g0358800, LOC_Os06g25250
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69KJ0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MESPPHATAS ADEMPSIWKE QHAQDAPPGF VPPMGPGEVA AVESLLGYEF RDKALVEEAL THGSFYYPYR PGVTYERLEY LGDAVLTCVV SREVFLTYGQ LQPGPLTRLR AANVDKEKLA RVAVVHGLHH FLRHKAPNLD GQITDFIEEL SMYPIHSNGL LDPPKVLCDV VESLIGAIYC DSNFNQEIWV QVFQKLADPL ISLETLGKHP VSELFECQK TRRGVKIVKD EWDKNLTVEV LIDGEMVGRA TYAQKKEIAQ NRAAKAALDK LKETLGQSQT EPMSAEVSEQ FNKIDLTGS
Source	Baculovirus
Target Names	Os06g0358800
Protein Names	Recommended name: Ribonuclease 3-like protein 3 EC= 3.1.26.- Alternative name(s): Ribonuclease III-like protein 3 Short name= RNase III-like protein 3
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.