



Recombinant *Oryza sativa* subsp. *japonica* BURP domain-containing protein 15 (BURP15)

Product Code	CSB-BP724354OFG
Abbreviation	BURP15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69QR8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	GR QMTAFAARTS PAAAAEAFWR AAMPGAPMPD AIVELLHHEH GVASAGGKAN GGGDGPPPPM NFNYDDYRAL PRSDAPSPDA LNRVAAVQNA DENGVSPPP PPPTVFFLED AVRVGESLPL PRPAADATAA GAAAATALPP LRLYTVRSVR AIEGSSFVVC RGETTAGAGV YGCRDAATGP ARAYAVDAAG GGGGDAVIAA VVCHADTSRW DPDHAAFRL GVRPGGAAVC RAVADAHILP TNKD
Source	Baculovirus
Target Names	BURP15
Protein Names	Recommended name: BURP domain-containing protein 15 Short name= OsBURP15 Alternative name(s): Anther-specific protein 8
Expression Region	19-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.