



Recombinant *Oryza sativa* subsp. japonica Homeobox-leucine zipper protein ROC8 (ROC8), partial

Product Code	CSB-MP724356OFG
Abbreviation	ROC8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69T58
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	ROC8
Protein Names	Recommended name: Homeobox-leucine zipper protein ROC8 Alternative name(s): GLABRA 2-like homeobox protein 8 HD-ZIP protein ROC8 Homeodomain transcription factor ROC8 Protein RICE OUTERMOST CELL-SPECIFIC 8
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.