



Recombinant *Oryza sativa* subsp. *japonica* Ribose-phosphate pyrophosphokinase 2, chloroplastic (Os06g0617800, LOC_Os06g41360)

Product Code	CSB-BP724360OFG
Abbreviation	Os06g0617800, LOC_Os06g41360
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69XQ6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	CVMGSEQV RVVVEEEGKT KKRMGVAEPR SAPPVWTPR APAQEARLAA LRTDGRDSRL KIFSGTANRP LAQEIASYLG VDLGKVLIKR FADGEIYVQL QESVRGCDVF LVQPTCSPVN ENLMELFVMI DACRRASARS ITVVIPYFGY ARADRKAQGR EAITAKLSAN LLTEAGSDRV IVCDIHSTQA LGYFDIPVDH IHGQPVILDY LASKTISKDL VVVSPDVGGV VRARAFAKKL SDAPLAIVDK RRQGHNMSEV MHLIGDVKGK VAIMVDDMID TAGTITSAAA LLKQEGAEAV YACSTHAVFS PPAIERLSGG IFEEVIVTNS ILLPEHKCFP QLTVLSMANL VAETIWHVHR DGSVSSIFQ
Source	Baculovirus
Target Names	Os06g0617800
Protein Names	Recommended name: Ribose-phosphate pyrophosphokinase 2, chloroplastic EC= 2.7.6.1 Alternative name(s): Phosphoribosyl pyrophosphate synthase 2
Expression Region	33-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.