



Recombinant Leifsonia xyli subsp. xyli Ribosomal RNA small subunit methyltransferase G (rsmG)

Product Code	CSB-BP724410LAAC
Abbreviation	rsmG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ABV7
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MTDSVELEPA AAAQLFGDRL PIARTFTANL ASQGEERGLI GPLELPRLWS RHILNSAIVA PLLRPGVVG D VGSGAGLPGL VLAIARPDVS FVLIEPMDRR VIWLNEQVAE LGLLNVDIVR ARAEDAKLAK PLDQVTARAV SAFRKLPLT APLLRDGGEL VLMKGAAAQA EIDGAAKEIS KFKVRGAEVI VLGEGVLDEV TRVIRATVR
Source	Baculovirus
Target Names	rsmG
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase G EC=2.1.1.- Alternative name(s): 16S rRNA 7-methylguanosine methyltransferase Short name= 16S rRNA m7G methyltransferase
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.