



Recombinant *Leifsonia xyli* subsp. *xyli* Ribosome-binding factor A (rbfA)

Product Code	CSB-EP724439LAAC
Abbreviation	rbfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AG48
Product Type	Recombinant Protein
Immunogen Species	<i>Leifsonia xyli</i> subsp. <i>xyli</i> (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MADPARARKL ADRIKVIVAE RLDRGLRDPR LGFVTITDVQ VTGDLQHATV FYTVYGTQEQE RADTAAALTA ATGMLRSEVG KNITARLTPT LQFQLDAIPE NAAHIEDLLR QARDHDTEVA GLAEGAAYAG EQDPYVKPRE AN
Source	E.coli
Target Names	rbfA
Protein Names	Recommended name: Ribosome-binding factor A
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.