



Recombinant Leifsonia xyli subsp. xyli Acetylglutamate kinase (argB)

Product Code	CSB-BP724441LAAC
Abbreviation	argB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AGD1
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MTTEDTFETA ERAAAIAKAA ALIESLPWLK TFHDRIIVVK FGGNAMVSEE LQRTFAEDMV YLRYAGLLPV IVHGGGPQIS AMLDRLGIQS EFRGGYRVT PEAMEVVRMV LTGQINRDIV AAINKHGPLA AGLSGEDAAL FQGRKRGAVV DGEAVDLGLV GDVIGVNPEA VLAQIDAGRI PVVSSIAPDI DEPGQSLNVN ADAAAAALAV ALRAAKLVIL TDVAGLYSDW PNRDSLISKI TARELRELLP GLES GMIPKM AACLA AVDGG VEKAEIIDGR IEHSILLEVF TQSGIGTEVA PV
Source	Baculovirus
Target Names	argB
Protein Names	Recommended name: Acetylglutamate kinase EC= 2.7.2.8 Alternative name(s): N-acetyl-L-glutamate 5-phosphotransferase NAG kinase Short name= AGK
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.