



Recombinant *Oryza sativa* subsp. *japonica* Splicing factor U2af small subunit B (U2AF35B)

Product Code	CSB-MP724487OFG
Abbreviation	U2AF35B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AUG0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAEHLASIFG TEKDRVNCPF YFKIGACRHG DRCSRLHNRP TVSPTIVLAN MYQRPDMITP GVDAQGQPID PEKMQEHFED FYEDIYEELS KFGEVETLNV CDNLADHMIG NVYVQFREEE QAVAAHNALQ GRFYSGRPII VEYSPVTDFR EATCRQFEEN SCNRGGYCNF MHVKQIGREL RRKLYGGRSR RSHGRSRSPS PRHRRGNRDR DDFRRERDGY RGGGDGYRGG GGGGGGDGYR GGDSYRGGGG GGRRGGGSRY DRYDDGGRRR HGSPRRRARS PVRESSEERR AKIEQWNRER EEKP
Source	Mammalian cell
Target Names	U2AF35B
Protein Names	Recommended name: Splicing factor U2af small subunit B Alternative name(s): U2 auxiliary factor 35 kDa subunit B U2 small nuclear ribonucleoprotein auxiliary factor small subunit B Short name= U2 snRNP auxiliary factor small subunit B
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.