



# Recombinant Gracilaria tenuistipitata var. liui Pyruvate dehydrogenase E1 component subunit beta (pdhB)

<b>Product Code</b>	CSB-BP724540GBAK
<b>Abbreviation</b>	pdhB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6B8T1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gracilaria tenuistipitata var. liui (Red alga)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTEVLMFDAL REATDEEMQN DSSVFILGED VGHYGGSYKV TKDLHSKYGD LRVLDTPIAE NSFMGMAIGA AITGLRPIVE GMNMSFLLLA FNQISNNAGM LRYTSGGNFQ IPIVIRGPGG VGRQLGAEHS QRLEAYFQAI PGLKIVACST PYNAGLLKS AIRDNNPVIF FEHVLLYNLK DELPNDEYFL PLDKAELVRD GLDVTILTYR RMRHHVMQAV VDLVNDGYNP EVIDLISLKP LDITSIAQSL MKTHKLIIVE ECMKTGGIGA EIIAQINDNY FDFLDAPIVR LSSQDIPTPY NGKLEKATVI YPQQIIEAVK SIV
<b>Source</b>	Baculovirus
<b>Target Names</b>	pdhB
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit beta EC= 1.2.4.1
<b>Expression Region</b>	1-323
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.