



# Recombinant Uncharacterized protein R06F6.12 (R06F6.12)

<b>Product Code</b>	CSB-MP724555CXY
<b>Abbreviation</b>	R06F6.12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6BEV5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDSEKENPAS PTTGREEEIP GSSPEGVYPA DEDHIFHEDQ APLRVESAKH EEEIVEQQQQ QPEDLEQGDM IVEDGDQQFM NMVQITQEDM YEAGFDVEAG FELNHLTEEQ LNIVVAISQQ RQAKQNEQEH NEEVVIEDGS HHHMVHHEMI ENQFEGDGNG PEEYDGNQ IIDNAMHIIL TNDGGVNITS KQKQFYVSPS EIANLNIDLN NLSTENVHQL VQLALPPIKE KAQDSAYNDQ APSTSYHHHH HEQLEAGKST RSPIIGETVQ IRTADGRLQD AVVKYVRGDS EYKIQLMNGE FAYATIDQML VPQRDRSDHE YQQVAPPVLL RRTDMASVRN AAQKRSANDD LCPPVLKKSQ QLAPVVDGPH LVHTPNFCCP ICDKKVYQKE PSYIVIRLPA CDSCTREKII VLDEQSS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	R06F6.12
<b>Protein Names</b>	Recommended name: Uncharacterized protein R06F6.12
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.