



Recombinant *Debaryomyces hansenii* Pre-mRNA-splicing factor SLU7 (SLU7)

Product Code	CSB-EP724582DIS
Abbreviation	SLU7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BMK7
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSSQEKHKPE YDSSGKEVNP YIPKFISAVP WYHNKSNDEK SDDYLSHQRS NVADEAKDHS IPQPGSGIND EFEIKGETEI KKVEDYDSKR DRWHGYEAQE WDKIAENWDK IKKKKQKTKN ASVEEDSDDT DYELELVELG LDSKDIKNNL KEDPLEKTIR DRQDVPAYIL NITSSNKIHY DPKSRLTKDP SKGFINDKNQ FVKKLTGEAK RLDNLQKFAW EQNRQQEEMK QREAFEQKLT GKKHSEGGPD EYQVDLNLNM EANPTAMMLQ ARHKQEQQQA SHDQKKS DLV AKYGGGEFLN KSKEFVNVTE SIQTGDKPIN KDKNGLKRSI YPEDNYS MNH QSIWGSYYSG GQWGYCCCKQ TTRNSRCTLN
Source	E.coli
Target Names	SLU7
Protein Names	Recommended name: Pre-mRNA-splicing factor SLU7
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.