



# Recombinant Debaryomyces hansenii Protein PXR1 (PXR1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP724612DIS                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Abbreviation</b>      | PXR1                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q6BUE3                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | GLAGTKVKQ RFGLDPRNTN WSNDTSRFGH QYLERMGWKP GKGLGLVEHA<br>TTSHVKVSIK DDNLGLGSKL AKKQKTDEFD SGECSGLDVF QRILGRLNGK<br>EDQINKEIER KRKDNIINGK WGIQFIKGEV LQSTWDKEAK KLLDKGTTKK<br>RKIDSVDSSE ESTTSSDSKQ HKKKKIKKDR KEKEEKKTEK ENSEIKKKKK<br>EKKEKKEKKE KKDKNEKKEK KDKNEKKEKK DKNEEKEKKE KKEKKEKKDK<br>KDKKDKKDKK EKKEVKEVTR DSMLMPKEQL NQQIATRLSV RSKWIKQKRA<br>SVMDAKALNE IFMVTs                                |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | PXR1                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: Protein PXR1 Alternative name(s): PinX1-related protein 1                                                                                                                                                                                                                                                                                                                                 |
| <b>Expression Region</b> | 2-316                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |