



Recombinant *Kluyveromyces lactis* Methionine aminopeptidase 2 (MAP2)

Product Code	CSB-EP724739KBK-B
Abbreviation	MAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CP80
Product Type	Recombinant Protein
Immunogen Species	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (<i>Candida sphaerica</i>)
Purity	>85% (SDS-PAGE)
Sequence	MVADTDVKEI VVDALPEVSE PAAVDDSEVT EDATVQDKKK KKKKKKKKGN NMKNIALIYP DEKYPEGQWM EYHQDFNLKR VTDEERRYLE RSEEYEQRGN DMRKGAEIHR RVRESVRNKI KPGMTLTEIA NLVEDGTRKF TGTDANGDHV DRPKSQGIAF PTGLSLNHCA AHFTPNAGDK TVLKFEDVMK VDFGVHVNGY IIDSFTIAF DPQYDNLLAA VKDATNTGIK EAGIDVRLTD IGAEIQEVME SYEVEINGET HQVKPCRNLG GHNINPYSIH GGKSVPIVKN GDNTKMEENE HFAIETFGST GRGYVIQEGE CSHYAKKPGS HPTPSLSSAK NLLKVIDENF GTIPFCRRYL DRLGEDKHVY ALNTLVRQGI VEDYPPLNDI KGSYTAQFEH TLILPHKKE IVSRGDDY
Source	E.coli
Target Names	MAP2
Protein Names	Recommended name: Methionine aminopeptidase 2 Short name= MetAP 2 EC= 3.4.11.18 Alternative name(s): Peptidase M 2
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.