



# Recombinant *Erwinia carotovora* subsp. atroseptica (Dimethylallyl)adenosine tRNA methylthiotransferase MiaB (miaB)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP724860EAAB-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Abbreviation</b>      | miaB                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | Q6D7K6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) ( <i>Erwinia carotovora</i> subsp. atroseptica)                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MTKKLHIKTW GCQMNEYDSS KMADLLGSTH GYELTEIAEE ADVLLLNTCS<br>IREKAQEKVF HQLGRWKALK DLNPNLIIGV GGCVASQEGA HIRERAHYVD<br>VIFGPQTLHR LPEMINHVQG TRSPIVDISF PEIEKFDRLP EPRADGPTAF<br>VSIMEGCNKY CTFCVVPYTR GEEVSRPCDD VLFEIAQLAA QGVREVNLLG<br>QNVNAYRGET YDGEICSFAE LLRLVAAIDG IDRVRFTTSH PIEFTDDIIS<br>VYEDTPELVS FLHLPVQSGS DRVLTMMKRR HTALEYKAI RKLHNARPGI<br>LISSDFIIGF PGETQADFEQ TMKLIADVNF DMSYSFVYSA RPGTPAADMV<br>DDVLEEEKKQ RLYLLQDRIT KQAMRFSRLM LGTVQRILVE GTSRKSMEL<br>SGRTENNRVV NFEGTPDMIG KFVDVEIVDV YTNSLRGIVV RTEDQMDLRV<br>HESPSSVIAR TRKENEIGVG FYQP |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Names</b>      | miaB                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: (Dimethylallyl)adenosine tRNA methylthiotransferase<br>MiaB EC= 2.-.-.- Alternative name(s): tRNA-i(6)A37 methylthiotransferase                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Expression Region</b> | 1-474                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                        |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.