



Recombinant *Erwinia carotovora* subsp. atroseptica Trigger factor (tig)

Product Code	CSB-MP724864EAAB
Abbreviation	tig
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D828
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MQVSVETTQG LGRRVTITVA ADIIESAVKS ELVNAAKKVR IDGFRKGKVP MNVVAQRYGA SVRQDVLGDL MQRNFVDAIL KEKINPVGAP NYTPGEYAVG GDFTYTVFED VYPEVELKGL EAIEVEKPVV EVTDADVDTM LDTLRKQQAT WKETDRAATA EDRATIDFTG SIDGEVFEGG KASDFVLAMG QSRMIPGFED GIVGHKAGEE FTINVNFPED YHAENLKGKA AQFVIVLKKV EERELPELTE EFIKRFGVAD GSQEGRLRAEV RKNMERELKG AVRNRVKTQV LDGLINANDI EVPVALIDGE IDVLRQAAQ RFGGNEKQAL ELPRELFEEQ AKRRVVIGLL LGEVISSNEL KADEARVNVL IEEMASAYED PQEVIEFYSK NKELLNNMRN VALEEQAVET VLAKAKVVEK SVSFNELMNQ TTAA
Source	Mammalian cell
Target Names	tig
Protein Names	Recommended name: Trigger factor Short name= TF EC= 5.2.1.8 Alternative name(s): PPlase
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.