



Recombinant *Erwinia carotovora* subsp. atroseptica Phosphopantetheine adenylyltransferase (coaD)

Product Code	CSB-MP724883EAAB
Abbreviation	coaD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DAV2
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MTTKAIYPGT FDPLTNGHLD LLTRAARLFD HVVLAIAASP SKNTLFTLDE RVALANDATQ HLSNVEVIGF SDLMVHFAQQ QKANILVRGL RAVADFEYEL QLAKMNHHLMP PTLESVFLMP SEEWSFISSS LVKEVARHGG DVTHFLPDAI VRALLEKLH
Source	Mammalian cell
Target Names	coaD
Protein Names	Recommended name: Phosphopantetheine adenylyltransferase EC= 2.7.7.3 Alternative name(s): Dephospho-CoA pyrophosphorylase Pantetheine-phosphate adenylyltransferase Short name= PPAT
Expression Region	1-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.