



Recombinant Danio rerio 28S ribosomal protein S15, mitochondrial (mrps15)

Product Code	CSB-MP724925DIL
Abbreviation	mrps15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGL8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	KQEFPSQLED LKTSMLKHEY ASVPLAQSVD DVVKLLTLE FANHSEKLRL KEEQLIQKVQ RDENDRSSTE VKVAILTARI RNFQEHLHKH PKDKANKRWM LMTIDKRKKL LKFLRRTRYD SFEKVCKELG ITYTFPPEYY RRVTRRWLAK KAFCKNVFKE VQRLKAAERE KQKAEAAERK KSSSSTNPQE TAA
Source	Mammalian cell
Target Names	mrps15
Protein Names	Recommended name: 28S ribosomal protein S15, mitochondrial Short name= MRP-S15 Short name= S15mt
Expression Region	71-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.