



Recombinant Danio rerio Coiled-coil domain-containing protein 65 (ccdc65)

Product Code	CSB-EP724929DIL-B
Abbreviation	ccdc65
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHI2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MPKKAGKKAG GKLAGLTEEE RLLYMQQKAQ AEEEIAKRKE DMLTHFLKDK LQKEEKNTVL NLHKLRQQWR AVLTQTKIAE LRNDMSVLSQ TFERVLDYKD NIIRSLQVDL SERELQSELK RSSHLHNVDC LLEIHKSR LA QLEFNFKSSL VELSSEYNTREQILSQHQQ ECVDLDNVMF SMEKHYSDDL GEAKRDYQST RNQIKKRNFD DKQAVKEQMD GVVEKLWQDQ QQVLNHYNES TRDKIITDD LMNKDVQSAK EIDSLKKHIQ KLQDSISTFR GQLSSGQTDK TAEQLRSERD ELAQEVQHFR VQLNAAQAIR KKHITKLT VQ SNDATKKLQE IVARGERLLR QCEMCCKLET EHEKVLPHYK SSLSEEEQKQ EKAKAMESSN EKLTQLMHDY SPLAKFWQRY NKVELDCLCV KREKLLLLQE NERLRLFLKQ YLDEVSVSDE SFRQQKLLVV SSPALQDTAA TDRHQQKRHV VQEAACIVQK QL
Source	E.coli
Target Names	ccdc65
Protein Names	Recommended name: Coiled-coil domain-containing protein 65
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.