



Recombinant Human Adenosine deaminase-like protein (ADAL)

Product Code	CSB-MP724932HU
Abbreviation	ADAL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHV7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MIEAEEQQPC KTDFYSELPK VELHAHLNGS ISSHTMKKLI AQKPDLKIHD QMTVIDKGGK RTLEECFQMF QTIHQLTSSP EDILMVTKDVI KEFADDGVK YLELRSTPRR ENATGMTKKT YVESILEGIK QSKQENLDID VRYLIAVDRR GGPLVAKETV KLAEEFFLST EGTVLGLDLS GDPTVGQAKD FLEPLLEAKK AGLKLALHLS EIPNQQKETQ ILLDLLPDRI GHGTFLNSGE GGLDLVDFV RQHRIPLELC LTSNVKSQTV PSYDQHFFGF WYSIAHPSVI CTDDKGVFAT HLSQEYQLAA ETFNLTQSQV WDLSEYISNY IFASDSTRSE LRKKWNHLKP RVLHI
Source	Mammalian cell
Target Names	ADAL
Protein Names	Recommended name: Adenosine deaminase-like protein EC= 3.5.4.-
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.