



Recombinant *Xenopus tropicalis* UPF0364 protein C6orf211 homolog

Product Code	CSB-MP724941XBF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DJA3
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MDPPCSLSAR FEGSFAYLTI RDRLPQILTK VIDTVHRNKN KFFEDNGEEG VEAEKKALSF FSKLRNEIQT DKPVLPLTDN QSDTELWNQY LDYQKTLLNK GETPSWFKSP WLYVECYMYR RIQEGVLVSP PISEYDVFRE GKTESFFQSQ PAIALCTYL QELKNNVANL SENQKKEELS KLLQVCLWGN KCDLSISGGL DNSQKFSILS SLESFMPFIL VNDMESVWAV LSGSKNLESG KELMKRVDIV LDNAGFELIT DFVLADALLS LRLASEVHFH AKCMPWFVSD TTKHDFNWTI KHLQAANKHW MSKCGVNWKE NLQRSCWIYH EHLFWTL PHE FCMMAQTAPD LYSELQKSDL VIFKGDLYNR KLTGDRKWEF TVPFFQALTT FHPAPLCSIR TLKADVQVGL KPGVGGEQLSS NEPDWMISGK YGVVQLSTSV
Source	Mammalian cell
Target Names	arnt1
Protein Names	Recommended name: UPF0364 protein C6orf211 homolog
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.