



Recombinant Influenza A virus Nucleoprotein (NP)

Product Code	CSB-MP724956IGJ
Abbreviation	NP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DPE6
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/Silky Chicken/Hong Kong/YU100/2002 H5N1 genotype X3)
Purity	≥85% (SDS-PAGE)
Sequence	MASQGTKRSY EQMETGGERQ NATEIRASVG RMVSGIGRFY IQMCTELKLS DYEGRLIQNS ITIERMVLSA FDERRNRYLE GHPSAGKDPK KTGGPIYRRR DGKWWRELIL YDKEEIRRIW RQANNGEDAT AGLTHLMIWH SNLNDATYQR TRALVRTGMD PRMCSLMQGS TLPRRSGAAG AAIKGVGTMV MELIRMIKRG INDRNFWRGE NGRRTRIAYE RMCNILKGKF QTAAQRAMMD QVRESRNPNG AEIEDLIFLA RSALILRGSV AHKSCLPACV YGLAVASGYD FEREGYSLVG IDPFRLQNS QVFSLIRPNE NPAHKSQVLW MACHSAAFED LRVSSFIRGT RVVPRGQLST RGVQIASNEN MEAMDSNTLE LRSRYWAIRT RSGGNTNQQR ASAGQISVQP TFSVQRNLPF ERATIMAAFT GNTEGRTSDM RTEIIRMES ARPEDVSFQG RGVFELSDEK ATNPVPSFD MNNEGSYFFG DNAEEYDN
Source	Mammalian cell
Target Names	NP
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.