



Recombinant Influenza A virus Matrix protein 1 (M)

Product Code	CSB-BP724958IHQ
Abbreviation	M
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DPQ6
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/Teal/China/2978.1/2002 H5N1 genotype W)
Purity	≥85% (SDS-PAGE)
Sequence	LLTEVETYVL SIIPSGPLKA EIAQKLEDVF AGKNTDLEAL MEWLKTRPIL SPLTKGILGF VFTLTPSER GLQRRRFVQN ALNGNGDPNN MDRAVKLYKK LKREITFHGA KEVALSYSTG ALASCMGLIY NRMGTVTTEV AFGLVCATCE QIADSQHRSH RQMATITNPL IRHENRMVLA STTAKAMEQM AGSSEQAAEA MEVANQARQM VQAMRTIGTH PNSSAGLRDN LLENLQAYQK RMGVQMQRFK
Source	Baculovirus
Target Names	M
Protein Names	Recommended name: Matrix protein 1 Short name= M1
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.