



Recombinant *Nymphaea alba* DNA-directed RNA polymerase subunit alpha (rpoA)

Product Code	CSB-EP725014NES-B
Abbreviation	rpoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6EW20
Product Type	Recombinant Protein
Immunogen Species	<i>Nymphaea alba</i> (White water-lily) (<i>Castalia alba</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MVREEVPVST RTLQWKCVES RADSKRLYYG RFVLSPLMKG QADTIGIAMR RALLGELEGT CITRAKSDKV PHEYSTVVGI EESVHEILMN LKKMVFRSDL YGTLDASICV RGPRHVTAQD IISPPSVEIV DTTQHIAVLT EPVDLCIELK IERGRGYCTR TQNNYQDGSY PIDAVSMPVR NANHSIHSYG NGNEKQEILF LEIWTNGSLT PKEALYEASR NLIDLFIPFL HAEQDINRN MEDNLKRASV PFFAFDDGLD NIKREILKR IFIDQLELPP RTYNCLKRSN IHTLLDLSK SQEDLMRIEH FRVEDVKQIF DILQKGFTID LLKNSNQFES R
Source	E.coli
Target Names	rpoA
Protein Names	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= PEP EC= 2.7.7.6 Alternative name(s): Plastid-encoded RNA polymerase subunit alpha Short name= RNA polymerase subunit alpha
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.