



Recombinant Uncharacterized protein

pXO2-37/BXB0038/GBAA_pXO2_0038 (pXO2-37, BXB0038, GBAA_pXO2_0038)

Product Code	CSB-BP725027BQE
Abbreviation	pXO2-37, BXB0038, GBAA_pXO2_0038
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6F047
Product Type	Recombinant Protein
Immunogen Species	Bacillus anthracis
Purity	≥85% (SDS-PAGE)
Sequence	MVVLKLGKNY PIDRDSRGRY RHSVFMRGA NVEVKLLNHE ELLKGKITTI EEYYCVLDVE GTSGSYQVTV NFSEVKYIKH EEFPTVEERS PKYSKGDKKT SFVFKIGEKI GCVFKDGKGI KGTLLEDAY YLYVKSEKEN YYTIMKGALS YIAHVKHEPL LLTNDFYTEE MKLADYKKPT EYVFSVGDTI TVYFPSGKNV SGVVLDESKY WLLQTEKRQ ITIFKGSYSY FKHETYETKA YLYVENRKLK KQLRSGDTN
Source	Baculovirus
Target Names	pXO2-37
Protein Names	Recommended name: Uncharacterized protein pXO2-37/BXB0038/GBAA_pXO2_0038
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.