



# Recombinant Salmonella paratyphi A Rhamnulokinase (rhaB)

<b>Product Code</b>	CSB-YP725116SBM
<b>Abbreviation</b>	rhaB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5PKG3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTFRHCVAVD LGASSGRVML ARYDSKHRTL TLREIHRFVN CLQKTDGFDT WDIDSLEKDI RLGLKKVCNE GILIDSIGID TWGVDYVLLD KQGQRVGLPV SYRDNRTTGI MPQALVQIGK SEIYRRSGIQ FLPFNTIYQL RALTKQQPEL TAQVAHALLM PDYFSYRLTG EMNWEYTNAT TTQLVNINTD DWDDTLLAWT GAKKSWFGRP SHPGNVIGDW ICPQGNRIPV VAVASHDTAS AVIASPLANK HSAYLSSGTW SLMGFESKKP YTTDEALAAN ITNEGGAEGR YRVLKNIMGL WLLQRVLKER RITDLPALIA QTEALPACRF LINPNDDRFI NPDDMHAEIQ AACRETDQPV PVSDAELARC IFDSLALLYA DILHELANLR GEKFTQLHIV GGGCQNALLN QLCADACGIR VMAGPVEAST LGNIGIQLMT LDELNNIDDF RQVVSANYDL TTYIPNPDSE IARHVAQFQP KRQTKELCA
<b>Source</b>	Yeast
<b>Target Names</b>	rhaB
<b>Protein Names</b>	Recommended name: Rhamnulokinase EC= 2.7.1.5 Alternative name(s): Rhamnulose kinase
<b>Expression Region</b>	1-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.