



Recombinant Idiomarina loihiensis Ribosomal RNA large subunit methyltransferase N (rlmN)

Product Code	CSB-BP725220IAAA
Abbreviation	rlmN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5QYC0
Product Type	Recombinant Protein
Immunogen Species	Idiomarina loihiensis (strain ATCC BAA-735 / DSM 15497 / L2-TR)
Purity	≥85% (SDS-PAGE)
Sequence	MTNAIDKKVN LLDLNREGIK EFFREMGEKP FRAEQVMKWL YHFCVDDFDE MTNLNKALRE KKLQVAEIRA PEVREQQQSS DGTIKFAMTL FDGQDVETVW IPEGDRATLC VSSQVGCALC CTFCTGAQG FNRNLSVAEI IGQVWRVNQL LGAYGKTGIK PVTNVMMGM GEPLNLNNV VPAMEMLDD LGFGLSKRRV TLSTSGVPA LEKLREIDV MLAISLHAPD DELRNEIVPI NKKYNIEEFL ASSRRYVEQS KAQRKVTVEY VMLDHVNDST DQAHALAKTL KDTPSKINLI PFNPFPGSDY GRSSNSRIDR FAKVLMYGL TVMVRKTRGD DIDAACGQLV GDVIDRTKRI LKRQKQKRG EAIKVT
Source	Baculovirus
Target Names	rlmN
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase N EC=2.1.1.192 Alternative name(s): 23S rRNA (adenine(2503)-C(2))-methyltransferase 23S rRNA m2A2503 methyltransferase
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.