



# Recombinant *Idiomarina loihiensis* Ribosomal RNA large subunit methyltransferase N (rlmN)

<b>Product Code</b>	CSB-EP725220IAAA-B
<b>Abbreviation</b>	rlmN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5QYC0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Idiomarina loihiensis</i> (strain ATCC BAA-735 / DSM 15497 / L2-TR)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTNAIDKKVN LLDLNREGIK EFFREMGEKP FRAEQVMKWL YHFCVDDFDE MTNLNKALRE KLKQVAEIRA PEVREQQQSS DGTIKFAMTL FDGQDVETVW IPEGDRATLC VSSQVGCALC CTF CSTGAQG FNRNLSVAEI IGQVWRVNQL LGAYGKTGIK PVTNVMMGM GEPLLNNV VPAMEMLDD LGFGLSKRRV TLSTSGVPA LEKLRERIDV MLAISLHAPD DELRNEIVPI NKKYNIEEFL ASSRRYVEQS KAQRKVTVEY VMLDHSVNDST DQAHALAKTL KDTPSKINLI PFNPFPGSDY GRSSNSRIDR FAKVLM EYGL TVMVRKTRGD DIDAACGQLV GDVIDRTKRI LKRQQKQRGG EIAVKTT
<b>Source</b>	E.coli
<b>Target Names</b>	rlmN
<b>Protein Names</b>	Recommended name: Ribosomal RNA large subunit methyltransferase N EC=2.1.1.192 Alternative name(s): 23S rRNA (adenine(2503)-C(2))-methyltransferase 23S rRNA m2A2503 methyltransferase
<b>Expression Region</b>	1-378
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.