



Recombinant Idiomarina loihiensis Outer-membrane lipoprotein carrier protein (lolA)

Product Code	CSB-BP725233IAAA
Abbreviation	lolA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R0C0
Product Type	Recombinant Protein
Immunogen Species	Idiomarina loihiensis (strain ATCC BAA-735 / DSM 15497 / L2-TR)
Purity	≥85% (SDS-PAGE)
Sequence	SG DESDKQALQS LLKQTQSLSA EFEQQVKDEQ GEVLQTLSGT LKLKRPANLY WHTKEPDESV MVANGKKVWY YNPFVEQVTI YAQQDMVDDS PLLLVLNSNG NQWQNYNVSF RDNRYFVEHQ TNGSKLELRF TDEKLTEITM VQAQGERTEL MLNNVALNET ISDEQFVFDV PADVDVDDQS
Source	Baculovirus
Target Names	lolA
Protein Names	Recommended name: Outer-membrane lipoprotein carrier protein
Expression Region	19-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.