



Recombinant Pongo abelii Saccharopine dehydrogenase-like oxidoreductase (SCCPDH)

Product Code	CSB-BP725270PYX
Abbreviation	SCCPDH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R5C9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	ATEQRPFHL VVFGASGFTG QFVTEEVARE QVDPERSSRL PWAVAGRSRE KLQRVLEKAA LKLGRPTLSS EVGIIICDIA NPASLDEMAK QATVVPNCVG PYRFYGEPIV KACIENGASC IDISGEPQFL ELMQLKYHEK AADKGVYIIG SSGFDSIPAD LGVIYTRNKM NGTLTAVESF LTIHSGPEGL SIHDGTWKSA IYGFQDQSNL RKLNRASNLK PVPLVGPGLK RRWPISYCRE LKGYSIPFMG SDVSVVRRRTQ RYLYENLEES PVQYAAAYATV GGITSVIKLM FAGLFFLFFV RFGIGRQLLI KFPWFFSFGY FSKQGPTQKQ IDAASFLLTF FGQGYSSQIG TDMKPNIKI CTQVKGPEAG YVATPIAMVQ AAMTLLNDAS HLPKAGGVFT PGAAFSKTKL IDRLNKHGIE FSVISSSEV
Source	Baculovirus
Target Names	SCCPDH
Protein Names	Recommended name: Saccharopine dehydrogenase-like oxidoreductase EC=1.-.-.
Expression Region	2-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.