



Recombinant Pongo abelii Sorting nexin-6 (SNX6)

Product Code	CSB-YP725284PYX
Abbreviation	SNX6
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5R613
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	MMEGLDDGPD FLSEEDRGLK AINVDLQSDA ALQVDISDAL SERDKVKFTV HTKSSLPNFK QNEFSVVRQH EEFIWLHDSF VENEDYAGYI IPPAPPRPDF DASREKLQKL GEGEGSMTKE EFTKMKQELE AEYLAIFKKT VAMHEVFLCR VAAHPILRKD LNFHVFLFLEYN QDLSVRGKNK KEKLEDFFKN MVKSADGVIV SGVKDVDDFF EHERTFLLEY HNRVKDASAK SDRMTRSHKS AADDYNRIGS SLYALGTQDS TDICKFFLKV SELFDKTRKI EARVSAEDL KLSDLLKYLL RESQAAKDLL YRRSRSLVDY ENANKALDKA RAKNKDVLQA ETSQQLCCQK FEKISESAKQ ELIDFKTRRV AAFRKNLVEL AELELKHAKG NLQLLQNCLA VLNGDT
Source	Yeast
Target Names	SNX6
Protein Names	Recommended name: Sorting nexin-6
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>