



Recombinant Pongo abelii Ermin (ERMN)

Product Code	CSB-BP725292PYX
Abbreviation	ERMN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R6D6
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MTDVPATFTQ AECNGDKPPE NGQQPITKIS EELTDVDSPL PHYRVEPSLE GAPTKGSQEE RRKLQGNMLL NSSMEEKILK ENPEEKLFVV HKAITDLSLQ ETSADEMTR EGRQWEKIPL SGSNQEIRRQ KERITEQPLK EEEDEDKKNK GHQAAEIEWL GFRKPSQADM LSHKHDEEQK VWDEEIDDDD DDNCNDDDEDE VRVIEFKKKH EEVSQFKEEG DASEDSPLSS ASSQAVTPDE QPTLGKKSDI SRNAYSRYNT ISYRKIRKGN TKQRIDEFES MMHL
Source	Baculovirus
Target Names	ERMN
Protein Names	Recommended name: Ermin Alternative name(s): Juxtanodin Short name= JN
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.