



Recombinant Pongo abelii Probable RNA-binding protein 18 (RBM18)

Product Code	CSB-YP725299PYX
Abbreviation	RBM18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R6W3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MEAETKTLPL ENASILSEGS LQEGHRLWIG NLDPKITEYH LLKLLQKFGK VKQFDLFLHK SGALEGQPRG YCFVNFETKQ EAEQAIQCLN GKLALSCKLV VRWAHAQVKR YDHNKNDKIL PISLEPSSST EPTQSNLSVT AKIKAIEAKL KMMANPDAE YPAAPVYSYF KPPDKKRTP YSRTAWKSRR
Source	Yeast
Target Names	RBM18
Protein Names	Recommended name: Probable RNA-binding protein 18 Alternative name(s): RNA-binding motif protein 18
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.