



# Recombinant Pongo abelii Peptidyl-prolyl cis-trans isomerase CWC27 homolog (CWC27)

<b>Product Code</b>	CSB-EP725318PYX
<b>Abbreviation</b>	CWC27
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R7W3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SNIYIQEPP TNGKVLLKTT AGDIDIELWS KEAPKACRNF IQLCLEAYYD NTIFHRVVPF FIVQGGDPTG TGSGGESIYG APFKDEFHSR LRFNRRGLVA MANAGSHDNG SQFFFTLGRA DELNNKHTIF GKVTGDTVYN MLRLSEVDID DEERPHNPHK IKSCEVLFNP FDDIIPREIK RPKKEKPEEE VKKLKPKGTK NFSLLSFGEE AEEEEEEVNR VSQSMKGKSK SSHDLLKDDP HLSSVPVVES EKGDAAGDLD DGGEGESAEH DEYIDGDEKN LMRERIAKKL KKDTSANVKS TGEGEVEKKS VNRSEELRKE ARQLKRELLA AKQKKVENAA KPAEKRSEEE EATPDGAVAE YRREKQKYEA LRKQQSKKGT SREDQTLALL NQFKSKLTQA IAETPENDIP ETEVEDDEGW MSHVLQFEDK SRKVKDASMQ DSDTFEIIDP RNPVNKRRRE ESKKLMREKK ERR
<b>Source</b>	E.coli
<b>Target Names</b>	CWC27
<b>Protein Names</b>	Recommended name: Peptidyl-prolyl cis-trans isomerase CWC27 homolog Short name= PPlase CWC27 EC= 5.2.1.8
<b>Expression Region</b>	2-473
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.