



Recombinant Pongo abelii Protein CWC15 homolog (CWC15)

Product Code	CSB-MP725443PYX
Abbreviation	CWC15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RE65
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	TTAARPTFE PARGGRGKGE GDLSQLSKQY SSRDLPSHTK IKYRQTTQDA PEEVRNRDFR RELEERERAA AREKNRDRPT REHTTSSSVS KKPRLDQIPA ANLDADDPLT DEEDEDFFEE SDDDDTAALL AELEKIKKER AEEQARKEQE QKAEERIRM ENILSGNPLL NLTGPSQPQA NFKVKRRWDD DVVFKNCAKG VDDQKKDKRF VNDTLRSEFH KKFMEKYIK
Source	Mammalian cell
Target Names	CWC15
Protein Names	Recommended name: Protein CWC15 homolog
Expression Region	2-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.