



Recombinant Pongo abelii WD repeat-containing protein 13 (WDR13)

Product Code	CSB-MP725459PYX
Abbreviation	WDR13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RF24
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MAAVWQQVLA VDARYNAYRT PTFPQFRTQY IRRRSQLLRE NAKAGHPPAL RRQYLRLRGQ LLGQRYGPLS EPGSARAYSN SIVRSSRTTL DRMEDFEDDP RALGARGHRR SVSRGSYQLQ AQMNRAVYED RPPGSVVPTS AAEASRAMAG DTLSSENYAF AGMYHVFDQH VDEAVPRVRF ANDDRHLAC CSLDGSISLC QLVPAPPTVL RVLRGHTRGV SDFAWSLSND ILVSTSLDAT MRIWASEDGR CIREIPDPDG AELLCTFQP VNNNLTVVGN AKHNVHVMNI STGKKVKGGS SKLTGRVLAL SFDAPGRLLW AGDDRGSVFS FLFDMATGKL TKAKRLVVHE GSPVTSISAR SWVSREARDP SLLINACLNK LLLYRVVDNE GTLQLKRSFP IEQSSHPVRS IFCPLMSFRQ GACVVTGSED MCVHFFDVER AAKAAVNKLQ GHSAPVLDVS FNCDESLAS SDASGMVIVW RREQK
Source	Mammalian cell
Target Names	WDR13
Protein Names	Recommended name: WD repeat-containing protein 13
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.