



Recombinant Rat Protein FAM110C (Fam110c)

Product Code	CSB-MP725493RA
Abbreviation	Fam110c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RKJ0
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MRALPTLDSL ARTRPSLGDS RAAEGTLTRQ PANKSAVERL AADRAKYVRS TLGSSRGPVS ENRVPEAPGV QHRNPIPSTL APAPVARRAV ARKPLRPDSL VIYRQKCEFV RGSADSSRV GLMKKFFQGS GKDKMSVTPE TTRVADEYKT KETEATWTKA SQAAAAPAS MLPPPAPVVA MESPALPFEV APRVPVRRSG VELRASRSKG LQRSQSDLSS RYSTARAESD TFFQYCGLDP EVVESLGREN FSAASDCVTP KMRSVSMATS DSSFSRHSED GLQEEELLEQ VPSTTSVVER NARIIKWLF CKKAKETPSQ KLQEPA
Source	Mammalian cell
Target Names	Fam110c
Protein Names	Recommended name: Protein FAM110C
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.