



Recombinant Pig Interleukin-27 subunit alpha (IL27)

Product Code	CSB-BP725497PI
Abbreviation	IL27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5S1V9
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	≥85% (SDS-PAGE)
Sequence	FP RPPGRPPLSL QELQREFKVS LHLARKLLSE VRVQARHFAE SHLPGVNLDL LPLGEQLPNV SLNFQAWRGL SDPERLCFLS MTLRPFHTLL GGLGSQGFWT SSERMQLWAI RDLRLDLQQH LRFQVLAAGF NLPGQEEEEEN EAGRELLPGA PGGPSKPAAQ VSWPRLLYTY QLLHSLELVL SRAMRDFLLL SRAGNPAPAL GFPTPSSPP
Source	Baculovirus
Target Names	IL27
Protein Names	Recommended name: Interleukin-27 subunit alpha Short name= IL-27 subunit alpha Short name= IL-27-A Short name= IL27-A Alternative name(s): p28
Expression Region	29-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is one of the subunits of a heterodimeric cytokine complex. This protein is related to interleukin 12A (IL12A). It interacts with Epstein-Barr virus induced gene 3 (EBI3), a protein similar to interleukin 12B (IL12B), and forms a complex that has been shown to drive rapid expansion of naive but not memory CD4(+) T cells. The complex is also found to synergize strongly with interleukin 12 to trigger interferon gamma (IFNG) production of naive CD4(+) T cells. The biological effects of this cytokine are mediated by the class I cytokine receptor (WSX1/TCRR).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.