



# Recombinant Pig Interleukin-27 subunit alpha (IL27)

<b>Product Code</b>	CSB-MP725497PI
<b>Abbreviation</b>	IL27
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5S1V9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	FP RPPGRPPLSL QELQREFKVS LHLARKLLSE VRVQARHFAE SHLPGVNLDL LPLGEQLPNV SLNFQAWRGL SDPERLCFLS MTLRPFHTLL GGLGSQGFWT SSERMQLWAI RLDLRDLQQH LRFQVLAAGF NLPGQEEEEEN EAGRELLPGA PGGPSKPAAQ VSWPRLLYTY QLLHSLELVL SRAMRDFLLL SRAGNPAPAL GFPTPSSPP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IL27
<b>Protein Names</b>	Recommended name: Interleukin-27 subunit alpha Short name= IL-27 subunit alpha Short name= IL-27-A Short name= IL27-A Alternative name(s): p28
<b>Expression Region</b>	29-239
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is one of the subunits of a heterodimeric cytokine complex. This protein is related to interleukin 12A (IL12A). It interacts with Epstein-Barr virus induced gene 3 (EBI3), a protein similar to interleukin 12B (IL12B), and forms a complex that has been shown to drive rapid expansion of naive but not memory CD4(+) T cells. The complex is also found to synergize strongly with interleukin 12 to trigger interferon gamma (IFNG) production of naive CD4(+) T cells. The biological effects of this cytokine are mediated by the class I cytokine receptor (WSX1/TCRR).
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.