



# Recombinant *Oryza sativa* subsp. japonica

## Probable protein phosphatase 2C 54

### (Os06g0179700, LOC\_Os06g08140)

<b>Product Code</b>	CSB-EP725538OFG
<b>Abbreviation</b>	Os06g0179700, LOC_Os06g08140
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5SMK6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCVEESEGAERLDFGEPAAA AADAGKSKSK SPDELPSPRM ERVCENTTAA DFKQNKSGNF VPNIRSGDWS DIGGRQYMED THVCITDLAK NFGYQSVDNE AISFYGVFDG HGGKDAAHFV RDNLPRIIVE DADFPLELEK VVRRSFVHAD NQFAKTTLSS GTTALTAMIF GRTLLIANAG DCRAVLSRCG TAIEMSVDHR PCSLSEKLRV ESLGGYVDDG YLNGLLGVTR ALGDWHLEGM KEAGNPGGGL SAPELKMITE LTKDDEFLLI GSDGIWDVFS NQNVVDFARR RLQEHNDVKS CCREIVEEAI KRGATDNLTA VLVSFHLEAP PQVRVSRPGR VARSISAEGL NSLRLLRNQ
<b>Source</b>	E.coli
<b>Target Names</b>	Os06g0179700
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 54 Short name= OsPP2C54 EC= 3.1.3.16
<b>Expression Region</b>	1-360
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.