



Recombinant Mouse Coiled-coil domain-containing protein 85A (Ccdc85a)

Product Code	CSB-MP725540MO
Abbreviation	Ccdc85a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SP85
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSKAAGGSAP AAESCPSAPA GASTPTGVDD LSKVTDEELL QWSKEELIRS LRRAEAEKVS AMLDHSNLIR EVNRRLQLHL GEIRGLKDIN QKLQEDNQEL RDLCCFLDDD RQKGKRVSR WQRLGRYTAG VMHKEVALYL QKLKELEVKQ EEVVKENMEL KELCMLLDEE KGVGCAGSRC SIDSQASLCQ LVSATPYVR DVGDGSSTSS TGSTDSPDHH KHHASGGSP EHLQKPRSEGS PEHTKHRSTS PEHLHKPRAS GTPDHSKALK GPSPEHHKPL CKGSPEQQRH PHPGSSPEVL PKHVLGSPE HFQKHRPGGS PEHTRHSGGS PEHLQKHALG GSLEHLPRAR GTSPEHLKQH YGASPDHKHA SGSGGSGGGS REGTLRRPAQ EDSSSHHRNV YSGMNESTLS YVRQLEARVR QLEENRMLP QGSFRLSSGA DGNSSLNSP ASFSGHTTSP QQPEPVVHSL KVVWRKLGDA AGSCPGIRQH LSGNQYKGPM
Source	Mammalian cell
Target Names	Ccdc85a
Protein Names	Recommended name: Coiled-coil domain-containing protein 85A
Expression Region	1-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.