



# Recombinant Mouse Coiled-coil domain-containing protein 85A (Ccdc85a)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP725540MO                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Abbreviation</b>      | Ccdc85a                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                               |
| <b>Uniprot No.</b>       | Q5SP85                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Sequence</b>          | MSKAAGGSAP AAESCPSAPA GASTPTGVDD LSKVTDEELL QWSKEELIRS<br>LRRAEAEKVS AMLDHSNLIR EVNRRLQLHL GEIRGLKDIN QKLQEDNQEL<br>RDLCCFLDDD RQKGKRVSR WQRLGRYTAG VMHKEVALYL QKLKELEV KQ<br>EEVVKENMEL KELCMLLDEE KGVGCAGSRC SIDSQASLCQ LVSATPYVR<br>DVGDGSSTSS TGSTDSPDHH KHHASGGSP E HLQKPRSEGS PEHTKHRSTS<br>PEHLHKPRAS GTPDHSKALK GPSPEHHKPL CKGSPEQQRH PHPGSSPEVL<br>PKHVLGSPE HFQKHRPGGS PEHTRHSGGS PEHLQKHALG GSLEHLPRAR<br>GTSPEHLKQH YGASPDHKHA SGSGGSGGGS REGTLRRPAQ<br>EDSSSHHRNV YSGMNESTLS YVRQLEARVR QLEENRMLP QGSFRLSSGA<br>DGNSSLNSP ASFSGHTTPS QQPEPVVHSL KVVWRKLGDA AGSCPGIRQH<br>LSGNQYKGPM |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Target Names</b>      | Ccdc85a                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: Coiled-coil domain-containing protein 85A                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Expression Region</b> | 1-500                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                      |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                                |



of lyophilized form is 12 months at -20°C/-80°C.