



# Recombinant Human Dapper homolog 2 (DACT2)

<b>Product Code</b>	CSB-EP725552HU
<b>Abbreviation</b>	DACT2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5SW24
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MWTPGGPPGS AGWDRRRLGA RLRAAFAGLQ ELQGLRATQQ ERVRGALALQ PPPAPAAPCG PHGLHGPEQQ LEAALAALQE QLSRLRQQDI GLKTHLDQLD LQISKLQLDV GTASGEALDS DSRPSSGFYE MSDGGSCSL TSCASVCS DH ISPSLGSLP VAQAHKARPS MGDWRPRSVD ETTVPAWRPQ ATEEGARPPG SVEDAGQPWG TFWPRPVSTG DLDRALPADT GLQKASADAE LLGLLCQGV D IPLHVPDPKY RQDLVSQGGR EVYPYPSPLH AVALQSPLFV LTKETPQRGG PSFPRESPRG PAGLNTIQTG PVLEAGPARA RAYIDRLLHL WGRETPAKGS EGEQGPLRHA ASPSPQRQGG WSTDGGGRL V FAPGREDEG GPAQSRGAGR GGPQQQGYMP LEGPQQSGSL PEEGSKPSNS CVLRETMVQA SPSSKAQQTP SAQDYGRGNI ISPSRMLDKS PSPASGHFAH PSFAASLKMG PPKSKAEKIK RSPMDKVLRF ARQPLLLDR PEGAHAAPQP SLEWDPAHWP TGRGGLQRRP ALAWEAPGRS CSESTLYPMP VLVPLAVAPQ ESHRTSAQAL FPFEASLLTS VARRKHRRWQ STVEISARAR LASCPEENLG PPRPVARRAG GPLARGRPSL VRQDAYTRSD SEPSKHAEC DPRFPSVIPE TSEGESSDHT TNRFGDRESS SSDEEGGAQS RDCDLALGYV AAGHAELAWT QEAPVSSGPL LSPVPKLCRI KASKALKKKI RRFQPTALKV MTMV
<b>Source</b>	E.coli
<b>Target Names</b>	DACT2
<b>Protein Names</b>	Recommended name: Dapper homolog 2
<b>Expression Region</b>	1-774
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.