



Recombinant Human Polynucleotide 5'-hydroxyl-kinase NOL9 (NOL9)

Product Code	CSB-MP725557HU
Abbreviation	NOL9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SY16
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ADSGLLLKR GSCRSTWLRV RKARPQLILS RRPRLRLGSL RWCGRRLLRW RLLQAQASGV DWREGARQVS RAAAARRPNT ATPSPIPSPT PASEPESEPE LESASSCHRP LLIPPVRPVG PGRALLLLPV EQGFTFSGIC RVTCLYGQVQ VFGFTISQGG PAQDIFSVYT HSCLSIHALH YSQPEKSKKE LKREARNLLK SHLNLDDRRW SMQNFSPQCS IVLLEHLKTA TVNFITSYPG SSIYIFVQESP TPQIKPEYLA LRSVGIRREK KRKGLQLTES TLSALEELVN VSCEEVDGCP VILVCGSQDV GKSTFNRYLI NHLLNSLPCV DYLECDLGQT EFTPPGCISL LNITEPVLGP PFTHLRTPQK MVYYGKPSCK NNYENYIDIV KYVFSAYKRE SPLIVNTMGW VSDQGLLLLI DLIRLLSPSH VVQFRSDHSK YMPDLTPQYV DDMDGLYTKS KTKMRNRRFR LAAFADALEF ADEEKESPVE FTGHKLIGVY TDFAFRITPR NRESHNKILR DLSILSYLSQ LQPPMPKPLS PLHSLTPYQV PFNAVALRIT HSDVAPTHIL YAVNASWVGL CKIQDDVRGY TNGPILLAQT PICDCLGFGI CRGIDMEKRL YHILTPVPPE ELRTVNCLLV GAIAIPHCVL KCQRGIEGTV PYVTTDYNFK LPGASEKIGA REPEEAHKEK PYRRPKFCRK MK
Source	Mammalian cell
Target Names	NOL9
Protein Names	Recommended name: Polynucleotide 5'-hydroxyl-kinase NOL9 EC= 2.7.1.- Alternative name(s): Nucleolar protein 9
Expression Region	2-702
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.